RHINOSEDAL®

INDICATIONS

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RHINOSEDAL®

INDICATIONS

RHINOSEDAL is a medical device indicated for adults and children from 2 years old, to decongest the nose by an osmotic action, in case of colds, allergic or non-allergic rhinitis and sinusitis.

COMPOSITION

Natural sea water (hypertonic: 2,2% NaCl), tamarind extract, mallow extract, purified water, citric acid.

PROPERTIES OF RHINOSEDAL

RHINOSEDAL is a ready-to-use nasal spray which acts mechanically on the upper respiratory tract allowing:

- **Decongest the nose** by the osmotic action of hypertonic seawater. RHINOSEDAL helps clear the airways and evacuate infectious agents and allergens.
- **Protect against the risk of superinfections** by the film-forming and balsamic action of mallow and tamarind extracts.

HOW TO USE RHINOSEDAL

Remove the cap before use. During the first use, make 1 or 2 pressures to prime the pump. Tilt your head slightly to the side, insert the tip into one nostril, spray 2 times in each nostril. Wipe off any excess that may run from the nose with a clean tissue. Thoroughly clean the tip with hot water after use. Replace the cap securely on the mouthpiece.

RHINOSEDAL should be given to children with the help of an adult.

CONTRAINDICATIONS

Do not use in children under 2 years old.

Do not use in case of hypersensitivity or individual allergy to one or more of the components.

SIDE EFFECTS

If you get any side effects, talk to your doctor or pharmacist. This also applies to any side effects not mentioned in this leaflet. You can also report side effects to the manufacturer at the address provided in this leaflet and to the competent authority directly via the national reporting system: National Agency for the Safety of Medicines and Health Products (ANSM) - Website: [www .signalement-sante.gouv.fr](http://www.signalement-sante.gouv.fr)

PRECAUTIONS FOR USE

Keep out of the reach and sight of children.

If symptoms persist for more than 15 days, consult your doctor.

The vial should only be used by one person for hygienic reasons and to avoid the transmission of pathogens that could remain on the tip.

STORAGE AND DISPOSAL CONDITIONS

Store RHINOSEDAL at a temperature below 30 ° C.

See the expiration date indicated on the box. This expiration date refers to the intact product and stored correctly.

Do not throw away RHINOSEDAL via wastewater or household waste.

Ask your pharmacist to eliminate RHINOSEDAL. These measures will help protect the environment.

Leaflet publication date: 15/04/2021.

	Manufacturer : ELERTE, 181 Rue André Karman, 93300 Aubervilliers – France. France. Tél : +33 (0)1 48 34 75 03		Store at a temperature below 30°C
	Batch Number		Do not use if the pack is damaged
	Expiration Date		Special Instructions for use
	Manufacturing Date		Date of first affixing of the CE marking: 15/04/2021